

# Interface Wizard Backend API Documentation

**Version:** 2.0.0

**Date:** December 20, 2025

**Author:** Interface Wizard Development Team

**File:** `mainwithfastapi.py`

---

## Table of Contents

1. Executive Summary
  2. System Architecture
  3. Technology Stack
  4. Complete Architecture Diagrams
  5. API Endpoints Reference
  6. Code Structure Explained
  7. Data Flow & Processing
  8. Real-Time Features
  9. Integration Guide
  10. Testing Guide
  11. Troubleshooting
  12. Glossary
- 

## Executive Summary

What is Interface Wizard?

Interface Wizard is a healthcare integration platform that:

- Converts natural language into HL7 v2.x messages
- Processes CSV/Excel files containing patient data
- Generates valid HL7 messages automatically
- Sends messages to Mirth Connect via MLLP protocol
- Provides real-time progress updates during processing

## Key Features


- **AI-Powered Generation** - Uses OpenAI GPT-4 or fallback generator
- **Real-Time Progress** - Server-Sent Events (SSE) for live updates
- **Dashboard Statistics** - Track processed patients and success rates
- **Flexible CSV Parsing** - Handles various column naming conventions
- **Mirth Integration** - MLLP protocol for HL7 message transmission
- **Validation System** - 3-tier field validation (System, Contextual, Critical)
- **Dual Mode Operation** - API mode (FastAPI) + Console mode

## Who Should Use This Document?

- **Backend Developers** - Understanding the API implementation
  - **Frontend Developers** - Integrating with the API endpoints
  - **System Administrators** - Deploying and maintaining the system
  - **Healthcare IT** - Understanding HL7 message generation
  - **QA Engineers** - Testing the complete workflow
  - **Project Managers** - Understanding system capabilities
- 

## System Architecture

## High-Level Overview



## Component Responsibilities

Component	Responsibility	Technology
**Frontend**	User interface, file upload, real-time display	Angular 17 / React 18

Component	Responsibility	Technology
**FastAPI Backend**	REST API, HL7 generation, validation	Python 3.9+, FastAPI
**OpenAI API**	AI-powered HL7 message generation	GPT-4o-mini
**Mirth Connect**	HL7 message routing and transformation	Mirth Connect 4.x
**OpenEMR**	Patient data storage	MySQL 8.0+

---

## Technology Stack

### Backend Technologies

Language: Python 3.9+  
 Web Framework: FastAPI 0.109.0  
 ASGI Server: Uvicorn 0.27.0  
 Data Validation: Pydantic 2.5.3  
 HL7 Parsing: python-hl7 0.4.5  
 Data Processing: Pandas 2.1.4  
 Excel Support: openpyxl 3.1.2  
 AI Integration: OpenAI 1.10.0  
 File Upload: python-multipart 0.0.6

### Key Python Libraries Explained

#### \*\*FastAPI\*\*

- Modern, high-performance web framework
- Automatic API documentation (Swagger UI)
- Built-in validation with Pydantic
- Async/await support for real-time features

#### \*\*Uvicorn\*\*

- Lightning-fast ASGI server
- Handles concurrent connections
- WebSocket and SSE support

#### \*\*Pydantic\*\*

- Data validation using Python type hints
- Automatic JSON serialization/deserialization
- Clear error messages for invalid data

#### **\*\*python-hl7\*\***

- HL7 v2.x message parsing and validation
- Segment and field extraction
- Message structure validation

#### **\*\*Pandas\*\***

- CSV/Excel file reading
- Data transformation and cleaning
- Flexible column mapping


#### **\*\*OpenAI\*\***


- GPT-4o-mini integration
- Natural language to HL7 conversion
- Fallback to deterministic generator if unavailable

---


## Complete Architecture Diagrams

### 1. Real-Time CSV Upload Workflow





## 2. HL7 Message Generation Flow



- Tier 3: Critical Fields (User data - MUS)
  - ✓ MSH-9 (Message Type) = "ADT^A01"
  - ✓ PID-3 (Patient ID) = "PAT001"
  - ✓ PID-5 (Patient Name) = "Doe^John"
  - ✓ PID-7 (DOB) = "19900115"
  - ✓ PID-8 (Gender) = "M"
- Result: ■ VALIDATION PASSED
- STEP 5: Send to Mirth Connect (Optional)
- MLLP Envelope:  
<VT>HL7\_MESSAGE<FS><CR>
- Where:
  - <VT> = 0x0B (Start Block)
  - <FS> = 0x1C (End Block)
  - <CR> = 0x0D (Carriage Return)
- TCP Socket → localhost:6661  
■ Wait for ACK/NACK response
- Response: MSA|AA|abc-123-def-456
- Result: ■ Message accepted by Mirth

### 3. Server-Sent Events (SSE) Architecture

```
graph TD; A[FRONTEND: EventSource Connection] --> B[HTTP GET (Keep-Alive Connection)]; B --> C[BACKEND: SSE Endpoint]; C --> D[SSE Stream (text/event-stream)]; D --> E[NETWORK: SSE Events];
```

The diagram illustrates a sequence of events:

- FRONTEND:** EventSource Connection
- HTTP GET (Keep-Alive Connection)**
- BACKEND:** SSE Endpoint
- SSE Stream (text/event-stream)**
- NETWORK:** SSE Events

Each step is represented by a black square icon. Arrows indicate the flow from one step to the next.

```
■ (500ms pause)
■ data: {"current_step":4,"step_status":"Generated 1/10",...}
■ (500ms pause)
■ data: {"current_step":4,"step_status":"Generated 2/10",...}
■ ...
■ data: {"current_step":6,"status":"completed",...}
■
■
■
■
■
■ FRONTEND: EventSource.onmessage Handler
■
■ eventSource.onmessage = (event) => {
■   const data = JSON.parse(event.data);
■
■   // Update wizard UI
■   setActiveStep(data.current_step); // Highlight 1-6
■   setStatusMessage(data.step_status); // Show status
■   updateProgressBar(data.generated_count);
■
■   // Close connection when done
■   if (data.status === 'completed') {
■     eventSource.close();
■     showResults(upload_id);
■   }
■ };
■
```

## 4. In-Memory Data Storage Structure

```
■ PYTHON DICTIONARIES (In-Memory Storage)
■
■
■ upload_sessions: Dict[str, Dict[str, Any]]
■
■ {
■   "abc-123-def-456": {
■     "id": "abc-123-def-456",
■     "filename": "patients.csv",
■     "status": "processing", // or "completed", "error"
■     "current_step": 4, // 1-6
■     "step_status": "Generated 5/10 messages",
■     "total_patients": 10,
■     "created_at": "2025-12-20T14:00:00",
■     "completed_at": null, // or ISO timestamp
■     "trigger_event": "ADT-A01",
■     "send_to_mirth": true,
■
■     // Step 2: Parsed patient data
■     "parsed_patients": [
■       {
■         "row_number": 2,
■         "patient_id": "PAT001",
■         "first_name": "John",
■         "last_name": "Doe",
■         "patient_name": "John Doe",
■         "dob": "19900115",
■         "gender": "M",
■         "address": "123 Main St, Boston, MA 02101",
■         "selected": true
■       },
■       ...
■     ],
■
■     // Step 3: Selected patient IDs
■     "selected_patient_ids": ["PAT001", "PAT002", ...],
■   }
■ }
```

```

// Step 4: Generated HL7 messages
"generated_messages": [
  {
    "row_number": 2,
    "patient_id": "PAT001",
    "patient_name": "John Doe",
    "hl7_message": "MSH|^~\&|SMART_APP|...",
    "validation": {
      "is_valid": true,
      "missing_fields": [],
      "message": "Validation passed"
    },
    "status": "success",
    "mirth_sent": true,           // if send_to_mirth enabled
    "mirth_ack": "MSA|AA|..."   // ACK from Mirth
  },
  ...
],
// Step 5: Mirth send statistics
"mirth_successful": 9,
"mirth_failed": 1
},
"xyz-789-ghi-012": { ... }, // Another session
...
}

```

```

dashboard_stats: Dict[str, int]

{
  "total_processed": 150,          // Total patients
  "hl7_messages_generated": 148,  // Total messages created
  "successful_sends": 145,        // Successfully sent
  "failed_sends": 3              // Failed to send
}

// Updated automatically as processing completes
// Used by GET /api/dashboard/stats

```

Note: Data persists only while server is running  
 Restart = all sessions lost (acceptable for dev)  
 Production: Consider adding database persistence

## API Endpoints Reference

### Complete Endpoint List

Method	Endpoint	Tag	Purpose
**GET**	`/`	General	API information
**GET**	`/health`	General	Health check

Method	Endpoint	Tag	Purpose
**GET**	`/api/dashboard/stats`	Dashboard	Dashboard statistics
**GET**	`/api/dashboard/system-status`	Dashboard	System health
**POST**	`/api/upload`	Real-Time Upload	Upload CSV with real-time progress
**GET**	`/api/upload/{id}/stream`	Real-Time Upload	SSE progress stream
**GET**	`/api/upload/{id}/status`	Real-Time Upload	Current status (polling)
**GET**	`/api/upload/{id}/results`	Real-Time Upload	Final results
**POST**	`/api/generate-hl7`	HL7 Generation	Generate from text
**POST**	`/api/validate-hl7`	HL7 Validation	Validate HL7 message
**POST**	`/api/send-to-mirth`	Mirth Integration	Send single message
**POST**	`/api/upload-excel`	Bulk Processing	Legacy bulk upload
**GET**	`/api/supported-events`	Reference	Supported trigger events

## 1. GET /api/dashboard/stats

**Purpose:** Get dashboard statistics for UI tiles

**Request:**

```
GET /api/dashboard/stats HTTP/1.1
Host: localhost:8000
```

**Response:**

```
{
  "total_processed": 150,
  "hl7_messages": 148,
  "successful_sends": 145,
  "failed_sends": 3,
  "success_rate": 98.0
}
```

**Response Fields:**

- `total_processed` (integer) - Total patients processed across all sessions
- `hl7_messages` (integer) - Total HL7 messages generated
- `successful_sends` (integer) - Successfully sent to Mirth

- `failed_sends` (integer) - Failed to send to Mirth
- `success_rate` (float) - Success percentage (0-100)

### Use Case:

```
// Frontend: Poll every 5 seconds to update dashboard
setInterval(async () => {
  const stats = await fetch('/api/dashboard/stats').then(r => r.json());

  document.getElementById('total-tile').textContent = stats.total_processed;
  document.getElementById('messages-tile').textContent = stats.hl7_messages;
  document.getElementById('success-rate-tile').textContent = stats.success_rate + '%';
}, 5000);
```

---

## 2. GET /api/dashboard/system-status

**Purpose:** Check system health for dashboard indicators

### Request:

```
GET /api/dashboard/system-status HTTP/1.1
Host: localhost:8000
```

### Response:

```
{
  "openemr_connection": {
    "status": "Active",
    "last_sync": "2025-12-20T14:30:00.123456"
  },
  "hl7_parser": {
    "status": "Running"
  },
  "message_queue": {
    "status": "Ready",
    "pending": 0
  }
}
```

### Response Fields:

- `openemr_connection.status` - "Active" if Mirth reachable, "Offline" otherwise
- `openemrconnection.lastsync` - Current timestamp (ISO format)
- `hl7_parser.status` - "Running" if OpenAI available, "Limited" if using fallback
- `message_queue.status` - Always "Ready" (for future queue implementation)

- `message_queue.pending` - Always 0 (for future queue implementation)

## How It Works:

```
# Tests Mirth connection by attempting socket connection
sock = socket.socket(socket.AF_INET, socket.SOCK_STREAM)
sock.settimeout(2)
result = sock.connect_ex((MIRTH_HOST, MIRTH_PORT))
sock.close()

if result == 0:
    mirth_status = "Active" # Connection successful
else:
    mirth_status = "Offline" # Connection failed
```

## Use Case:

```
// Frontend: Show green/red indicators
const status = await fetch('/api/dashboard/system-status').then(r => r.json());

if (status.openemr_connection.status === 'Active') {
    showGreenDot('mirth-indicator');
} else {
    showRedDot('mirth-indicator');
}
```

---

## 3. POST /api/upload

**Purpose:** Upload CSV file and start background processing with real-time updates

### Request:

```
POST /api/upload HTTP/1.1
Host: localhost:8000
Content-Type: multipart/form-data; boundary=----Boundary

----Boundary
Content-Disposition: form-data; name="file"; filename="patients.csv"
Content-Type: text/csv

Patient ID,First Name,Last Name,DOB,Gender,Address
PAT001,John,Doe,1990-01-15,M,123 Main St
----Boundary
Content-Disposition: form-data; name="trigger_event"

ADT-A01
----Boundary
Content-Disposition: form-data; name="send_to_mirth"

true
----Boundary--
```

## Request Parameters:

- `file` (file, required) - CSV or Excel file with patient data
- `trigger_event` (string, optional) - HL7 trigger event (default: "ADT-A01")
- `sendtomirth` (boolean, optional) - Send to Mirth after generation (default: false)

## Response:

```
{  
  "upload_id": "abc-123-def-456",  
  "filename": "patients.csv",  
  "status": "processing",  
  "message": "Upload started. Connect to /api/upload/abc-123-def-456/stream for real-time  
updates"  
}
```

## Response Fields:

- `upload_id` (string) - Unique session ID (UUID) - **Save this for next steps!**
- `filename` (string) - Original filename
- `status` (string) - Always "processing" on successful upload
- `message` (string) - Instructions for connecting to SSE stream

## What Happens Next:

1. Session created with unique `upload_id`
2. CSV file read and parsed
3. Background async task starts (`processcsvwith_progress`)
4. Processing happens in background through 6 steps
5. Frontend should connect to `/api/upload/{upload_id}/stream` for updates

## Use Case:

```
// Frontend: Upload file and get upload_id  
async function uploadCSV(file) {  
  const formData = new FormData();  
  formData.append('file', file);  
  formData.append('trigger_event', 'ADT-A01');  
  formData.append('send_to_mirth', true);  
  
  const response = await fetch('/api/upload', {  
    method: 'POST',  
    body: formData  
  });  
  
  const data = await response.json();  
  console.log('Upload ID:', data.upload_id);  
  
  // Now connect to SSE stream  
  connectToProgressStream(data.upload_id);  
}
```

---

## 4. GET /api/upload/{upload\_id}/stream

**Purpose:** Stream real-time progress updates using Server-Sent Events (SSE)

**Request:**

```
GET /api/upload/abc-123-def-456/stream HTTP/1.1
Host: localhost:8000
Accept: text/event-stream
```

**Response: (Server-Sent Events stream)**

```
data: {"current_step":1,"step_status":"File uploaded successfully","total_patients":0,
"status":"processing"}

data: {"current_step":2,"step_status":"Parsing CSV data...","total_patients":10,"status":
"processing"}

data: {"current_step":3,"step_status":"Selected 10 patients","total_patients":10,"status":
"processing"}

data: {"current_step":4,"step_status":"Generated 1/10 messages","total_patients":10,
"status":"processing","generated_count":1}

data: {"current_step":4,"step_status":"Generated 2/10 messages","total_patients":10,
"status":"processing","generated_count":2}

data: {"current_step":4,"step_status":"Generated 5/10 messages","total_patients":10,
"status":"processing","generated_count":5}

data: {"current_step":5,"step_status":"Sending to Mirth Connect...","total_patients":10,
"status":"processing","mirth_successful":3,"mirth_failed":0}

data: {"current_step":6,"step_status":"Processing complete!","total_patients":10,"status":
"completed","upload_id":"abc-123-def-456","successful":10,"failed":0}
```

**Event Data Fields:**

- `current_step` (integer) - Current step number (1-6)
- `step_status` (string) - Human-readable status message
- `total_patients` (integer) - Total patient count
- `status` (string) - "processing", "completed", or "error"
- `generated_count` (integer, optional) - Messages generated (step 4)
- `mirth_successful` (integer, optional) - Successful Mirth sends (step 5)
- `mirth_failed` (integer, optional) - Failed Mirth sends (step 5)

- `upload_id` (string, final event) - Session ID
- `successful` (integer, final event) - Successful message count
- `failed` (integer, final event) - Failed message count
- `error` (string, if error) - Error message

**Update Frequency:** Every 500ms (0.5 seconds)

### Use Case:

```
// Frontend: Connect to SSE stream and update UI
function connectToProgressStream(uploadId) {
  const eventSource = new EventSource(`/api/upload/${uploadId}/stream`);

  eventSource.onmessage = (event) => {
    const data = JSON.parse(event.data);

    // Update wizard step indicator
    highlightStep(data.current_step); // Highlight step 1-6

    // Update status message
    document.getElementById('status').textContent = data.step_status;

    // Update progress bar (if in step 4)
    if (data.generated_count && data.total_patients) {
      const percent = (data.generated_count / data.total_patients) * 100;
      document.getElementById('progress-bar').style.width = percent + '%';
      document.getElementById('progress-text').textContent =
        `${data.generated_count}/${data.total_patients}`;
    }

    // When complete, close connection and show results
    if (data.status === 'completed') {
      eventSource.close();
      loadResults(uploadId);
    }
  }

  // Handle errors
  if (data.status === 'error') {
    eventSource.close();
    showError(data.error);
  }
};

eventSource.onerror = () => {
  console.error('SSE connection error');
  eventSource.close();
};

}
```

### React Example:

```
import { useEffect, useState } from 'react';

function UploadProgress({ uploadId }) {
  const [progress, setProgress] = useState({});

  useEffect(() => {
    const eventSource = new EventSource(`/api/upload/${uploadId}/stream`);

    eventSource.onmessage = (event) => {
      const data = JSON.parse(event.data);
      setProgress(data);
    };
  });
}
```

```

        if (data.status === 'completed' || data.status === 'error') {
          eventSource.close();
        }
      );
      return () => eventSource.close();
    }, [uploadId]);
  return (
    <div>
      <h3>Step {progress.current_step}/6</h3>
      <p>{progress.step_status}</p>
      {progress.generated_count && (
        <progress
          value={progress.generated_count}
          max={progress.total_patients}
        />
      )}
    </div>
  );
}

```

---

## 5. GET /api/upload/{upload\_id}/status

**Purpose:** Get current upload status (polling alternative to SSE)

**Request:**

```
GET /api/upload/abc-123-def-456/status HTTP/1.1
Host: localhost:8000
```

**Response:**

```
{
  "upload_id": "abc-123-def-456",
  "filename": "patients.csv",
  "status": "processing",
  "current_step": 4,
  "step_status": "Generated 5/10 messages",
  "total_patients": 10,
  "created_at": "2025-12-20T14:00:00.123456",
  "completed_at": null
}
```

**Response Fields:**

- **upload\_id** (string) - Session ID
- **filename** (string) - Original filename
- **status** (string) - "processing", "completed", or "error"

- `current_step` (integer) - Current step (1-6)
- `step_status` (string) - Status message
- `total_patients` (integer) - Total patient count
- `created_at` (string) - Upload timestamp (ISO format)
- `completed_at` (string | null) - Completion timestamp or null if still processing

### **Use Case (Polling Alternative):**

```
// For environments that don't support SSE
async function pollStatus(uploadId) {
  const interval = setInterval(async () => {
    const status = await fetch(`/api/upload/${uploadId}/status`)
      .then(r => r.json());

    // Update UI
    updateStepIndicator(status.current_step);
    updateStatusMessage(status.step_status);

    // Stop polling when complete
    if (status.status === 'completed' || status.status === 'error') {
      clearInterval(interval);
      loadResults(uploadId);
    }
  }, 1000); // Poll every 1 second
}
```

---

## **6. GET /api/upload/{upload\_id}/results**

**Purpose:** Get complete results after processing is finished

### **Request:**

```
GET /api/upload/abc-123-def-456/results HTTP/1.1
Host: localhost:8000
```

### **Response:**

```
{
  "upload_id": "abc-123-def-456",
  "filename": "patients.csv",
  "status": "completed",
  "current_step": 6,
  "total_patients": 10,
  "successful": 9,
  "failed": 1,
  "mirth_successful": 9,
  "mirth_failed": 0,
  "created_at": "2025-12-20T14:00:00.123456",
  "completed_at": "2025-12-20T14:00:23.789012",
  "messages": [
    {
      "text": "All patients processed successfully. Total: 10, Successful: 9, Failed: 1."}
```

```

    "row_number": 2,
    "patient_id": "PAT001",
    "patient_name": "John Doe",
    "hl7_message": "MSH|^~\\|
&|SMART_APP|SMART_FAC|REC_APP|REC_FAC|20251220140005||ADT^A01|abc-123|P|2.5\
nEVN|A01|20251220140005\nPID|1||PAT001||Doe^John||19900115|M||123 Main St, Boston, MA
02101\nPV1|1|I",
    "validation": {
        "is_valid": true,
        "missing_fields": [],
        "message": "Validation passed"
    },
    "status": "success",
    "mirth_sent": true,
    "mirth_ack": "MSA|AA|abc-123"
},
{
    "row_number": 3,
    "patient_id": "PAT002",
    "patient_name": "Jane Smith",
    "hl7_message": "MSH|...",
    "validation": {
        "is_valid": false,
        "missing_fields": ["PID-7"],
        "message": "Missing critical fields"
    },
    "status": "validation_failed"
}
],
"error": null
}

```

## Response Fields:

- `upload_id` (string) - Session ID
- `filename` (string) - Original filename
- `status` (string) - "completed" or "error"
- `current_step` (integer) - Always 6 if completed
- `total_patients` (integer) - Total patient count
- `successful` (integer) - Successfully generated messages
- `failed` (integer) - Failed to generate
- `mirth_successful` (integer) - Successfully sent to Mirth
- `mirth_failed` (integer) - Failed to send to Mirth
- `created_at` (string) - Upload timestamp
- `completed_at` (string) - Completion timestamp
- `messages` (array) - Array of all message results
- `error` (string | null) - Error message if status is "error"

## Message Object Fields:

- `row_number` (integer) - Excel row number (starts at 2)
- `patient_id` (string) - Patient ID from CSV

- `patient_name` (string) - Full name
- `hl7_message` (string) - Generated HL7 message
- `validation` (object) - Validation result
- `status` (string) - "success", "validation\_failed", or "error"
- `mirth_sent` (boolean, optional) - If sent to Mirth
- `mirth_ack` (string, optional) - ACK from Mirth
- `error` (string, optional) - Error message if generation failed

## Use Case:

```
// Frontend: Display results in Step 6 completion screen
async function showResults(uploadId) {
  const results = await fetch(`/api/upload/${uploadId}/results`)
  .then(r => r.json());

  // Show summary
  document.getElementById('total').textContent = results.total_patients;
  document.getElementById('successful').textContent = results.successful;
  document.getElementById('failed').textContent = results.failed;
  document.getElementById('mirth-sent').textContent = results.mirth_successful;

  // Show message table
  const tbody = document.getElementById('results-table');
  results.messages.forEach(msg => {
    const row = tbody.insertRow();
    row.innerHTML =
      ` ${msg.patient_name}</td>       <td>${msg.patient_id}</td>       <td>${msg.status === 'success' ? '✓' : '✗'}</td>       <td>${msg.mirth_sent ? 'Yes' : 'No'}</td>       <td><button onclick="showHL7('${escape(msg.hl7_message})'}">View</button></td>`;   }); } |
```

## Code Structure Explained

### File Organization

```
main_with_fastapi.py (1,391 lines)
■
■■■ Imports (Lines 1-36)
■■■ Standard library (os, sys, re, socket, time, etc.)
■■■ Third-party (FastAPI, Pydantic, Pandas, etc.)
■■■ OpenAI SDK (optional)
■
■■■ Configuration (Lines 45-65)
■■■ OPENAI_API_KEY
■■■ MIRTH_HOST, MIRTH_PORT
■■■ FIELD_TIERS (validation)
■■■ TRIGGER_EVENT_MAPPING
```

```

■■■ Pydantic Models (Lines 67-96)
    ■■■ HL7GenerationRequest
    ■■■ HL7ValidationResponse
    ■■■ HL7GenerationResponse
    ■■■ MirthSendResponse
    ■■■ BatchProcessingResponse

■■■ FastAPI App Initialization (Lines 98-146)
    ■■■ FastAPI app with metadata
    ■■■ Global client_wrapper
    ■■■ upload_sessions dictionary
    ■■■ dashboard_stats dictionary

■■■ Helper Functions (Lines 148-172)
    ■■■ _normalize_colname() - Column name normalization
    ■■■ _find_column() - Flexible column matching

■■■ Client Wrapper Class (Lines 174-222)
    ■■■ __init__() - Initialize OpenAI client
    ■■■ generate() - Main generation method
    ■■■ generate_via_api() - OpenAI API call

■■■ HL7 Generation (Lines 224-293)
    ■■■ fallback_hl7_generator() - Local generator
    ■■■ generate_hl7_message() - Main entry point

■■■ HL7 Validation (Lines 295-337)
    ■■■ validate_hl7_structure() - Basic parsing
    ■■■ validate_required_fields_api() - 3-tier validation

■■■ Mirth Integration (Lines 339-375)
    ■■■ send_to_mirth() - MLLP protocol sender

■■■ Batch Processing (Lines 377-465)
    ■■■ process_excel_batch() - Synchronous batch processing

■■■ Async Background Processing (Lines 467-675)
    ■■■ process_csv_with_progress() - 6-step async workflow

■■■ API Endpoints (Lines 677-1180)
    ■■■ startup_event() - Initialize on startup
    ■■■ root() - API info
    ■■■ health_check() - Health endpoint
    ■■■ generate_hl7_from_command() - Text to HL7
    ■■■ validate_hl7_message() - Validate HL7
    ■■■ send_hl7_to_mirth() - Send single message
    ■■■ upload_and_process_excel() - Legacy bulk upload
    ■■■ get_supported_trigger_events() - Reference data
    ■■■ get_dashboard_stats() - Dashboard statistics ■ NEW
    ■■■ get_system_status() - System health ■ NEW
    ■■■ upload_csv_with_realtime_progress() - Real-time upload ■ NEW
    ■■■ stream_upload_progress() - SSE stream ■ NEW
    ■■■ get_upload_status() - Status polling ■ NEW
    ■■■ get_upload_results() - Final results ■ NEW

■■■ Console Mode Functions (Lines 1182-1364)
    ■■■ upload_excel_file() - File dialog
    ■■■ validate_required_fields() - Console validation
    ■■■ display_hl7_details() - Display segments
    ■■■ console_main() - Main console loop

■■■ Main Entry Point (Lines 1366-1391)
    ■■■ main() - Argument parsing and mode selection

```

## Key Functions Explained

## 1. `process\_csv\_with\_progress()` (Lines 467-675)

**Purpose:** Background async task that processes CSV through 6 steps

### Parameters:

- `upload_id` (str) - Session ID to update
- `df` (DataFrame) - Parsed CSV data
- `trigger_event` (str) - HL7 trigger event (default: "ADT-A01")
- `sendlomirth_flag` (bool) - Whether to send to Mirth (default: False)

### Flow:

```
async def process_csv_with_progress(...):
    session = upload_sessions[upload_id] # Get session reference

    try:
        # Step 1: File uploaded (0.5s)
        session["current_step"] = 1
        session["step_status"] = "File uploaded successfully"
        await asyncio.sleep(0.5)

        # Step 2: Parse CSV (0.3s + parsing time)
        # - Map column names flexibly
        # - Extract patient data
        # - Normalize dates, addresses
        # - Store in parsed_patients array

        # Step 3: Select patients (0.5s)
        # - Auto-select all patients
        # - Store selected_patient_ids

        # Step 4: Generate HL7 messages
        # - Loop through each patient
        # - Call generate_hl7_message()
        # - Validate generated message
        # - Update progress every iteration
        # - 1 second delay between messages (rate limiting)

        # Step 5: Send to Mirth (if enabled)
        # - Loop through generated messages
        # - Send via MLLP protocol
        # - Track success/failure counts
        # - 0.2 second delay between sends

        # Step 6: Complete
        # - Mark status as "completed"
        # - Set completion timestamp
        # - Update dashboard statistics

    except Exception as e:
        session["status"] = "error"
        session["error"] = str(e)
```

### Why Async?

- Non-blocking: Server can handle other requests while processing
- Real-time updates: Can update session state while running

- Rate limiting: Uses `await asyncio.sleep()` for delays
- Background execution: Started with `asyncio.create_task()`

## 2. `send\_to\_mirth()` (Lines 339-375)

**Purpose:** Send HL7 message to Mirth Connect via MLLP protocol

### MLLP Protocol:

MLLP (Minimal Lower Layer Protocol) Envelope:

`<VT> + HL7_MESSAGE + <FS> + <CR>`

Where:

`<VT>` = 0x0B (ASCII 11) - Start Block  
`<FS>` = 0x1C (ASCII 28) - End Block  
`<CR>` = 0x0D (ASCII 13) - Carriage Return

### Implementation:

```
def send_to_mirth(hl7_message: str, host: str = "localhost", port: int = 6661):
    sock = None
    try:
        # Build MLLP envelope
        START_BLOCK = b"\x0b"
        END_BLOCK = b"\x1c\x0d"
        hl7_bytes = hl7_message.replace("\n", "\r").encode("utf-8")
        mllp_message = START_BLOCK + hl7_bytes + END_BLOCK

        # Create TCP socket
        sock = socket.socket(socket.AF_INET, socket.SOCK_STREAM)
        sock.settimeout(15) # 15 second timeout
        sock.connect((host, port))

        # Send message
        sock.sendall(mllp_message)

        # Wait for ACK/NACK
        response = sock.recv(4096)
        ack_message = response.decode("utf-8", errors="ignore").strip("\x0b\x1c\x0d")

        # Check ACK code
        if "AA" in ack_message or "CA" in ack_message:
            return True, ack_message # Success
        else:
            return False, ack_message # Rejected

    except socket.timeout:
        return False, "Connection timeout!"
    except ConnectionRefusedError:
        return False, "Connection refused!"
    finally:
        if sock:
            sock.close() # Always close socket
```

### ACK Codes:

- **AA** - Application Accept (success)

- CA - Commit Accept (success)
- AE - Application Error (rejected)
- AR - Application Reject (rejected)

### 3. `generate\_hl7\_message()` (Lines 279-293)

**Purpose:** Generate HL7 message from patient data

**Flow:**

```
def generate_hl7_message(client: ClientWrapper, command: str) -> str:
    # Try OpenAI first
    if client.has_remote:
        try:
            return client.generate_via_api(prompt)
        except Exception:
            # Fallback to deterministic generator
            return fallback_hl7_generator(command)
    else:
        # No OpenAI, use fallback
        return fallback_hl7_generator(command)
```

#### OpenAI Prompt:

```
prompt = f"""
You are an HL7 v2.x message generator. Generate a valid HL7 message based on
the following command:

{command}

CRITICAL REQUIREMENTS:
- Use HL7 v2.x pipe-delimited format.
- Auto-populate MSH-11='P', MSH-12='2.5', EVN-2='{current_time}'
- Only include segments explicitly required or implied.
- For ADTs include MSH, EVN, PID, PV1 minimum.
- Output ONLY the HL7 message text (no explanation).
"""


# ...
```

#### Fallback Generator:

- Parses command text using regex
- Extracts patient fields (ID, name, DOB, gender, address)
- Builds HL7 message deterministically
- Always generates valid HL7 v2.5 format

## Data Flow & Processing

## Complete CSV Processing Flow





# Flexible CSV Column Mapping

**Problem:** CSV files have different column naming conventions

## Examples:

- "First Name" vs "FirstName" vs "first\_name" vs "Given Name"
  - "DOB" vs "Date of Birth" vs "birthdate"
  - "Gender" vs "Sex"

**Solution:** Flexible column mapping with normalization

```
# Define multiple possible names for each field
candidate_map = {
    'patient_first_name': [
        "Patient First Name",
        "First Name",
        "Firstname",
        "First Name"
    ]
}
```

```

        "given_name",
        "first_name"
    ],
    'dob': [
        "DOB",
        "Date of Birth",
        "birthdate",
        "date_of_birth"
    ],
    'gender': [
        "Gender",
        "Sex"
    ]
}

# Normalize column names for comparison
def _normalize_colname(name: str) -> str:
    s = name.strip().lower()
    s = re.sub(r"[^\w\s]", " ", s) # Remove punctuation
    s = re.sub(r"\s+", " ", s)      # Multiple spaces → single space
    s = s.replace(" ", "_")        # Spaces → underscores
    return s

# Find best matching column
for key, candidates in candidate_map.items():
    found = _find_column(df, candidates)
    mapping[key] = found

# Extract patient data using mapped columns
first_name = row.get(mapping.get('patient_first_name'))

```

**Result:** Works with ANY reasonable CSV column naming!

---

## Real-Time Features

### Server-Sent Events (SSE) Explained

#### What is SSE?

Server-Sent Events (SSE) is a standard for real-time server-to-client communication:

#### Characteristics:

- **One-way:** Server pushes updates to client
- **HTTP:** Uses standard HTTP (not WebSocket)
- **Text-based:** Sends text data
- **Auto-reconnect:** Browser automatically reconnects if connection drops

- **Simple:** Easier to implement than WebSockets

## SSE vs Polling vs WebSockets

Feature	SSE	Polling	WebSockets
**Direction**	Server → Client	Client → Server	Bidirectional
**Protocol**	HTTP	HTTP	Custom protocol
**Complexity**	Low	Very low	High
**Reconnection**	Automatic	Manual	Manual
**Server load**	Low	High	Low
**Browser support**	All modern browsers	All browsers	All modern browsers
**Use case**	Progress updates, live feeds	Simple status checks	Chat, gaming, collaboration

### When to use what:

- **SSE:** Progress updates, notifications, live data feeds (our use case)
- **Polling:** Simple status checks, older browser support required
- **WebSockets:** Real-time bidirectional communication (chat, gaming)

## SSE Implementation in Code

### Backend (FastAPI):

```

@app.get("/api/upload/{upload_id}/stream")
async def stream_upload_progress(upload_id: str):
    async def event_generator():
        session = upload_sessions[upload_id]

        # Stream events while processing
        while session["status"] == "processing":
            # Build event data
            event_data = {
                "current_step": session.get("current_step", 0),
                "step_status": session.get("step_status", ""),
                "status": session["status"]
            }

            # Send SSE event (CRITICAL FORMAT!)
            yield f"data: {json.dumps(event_data)}\n\n"
            #           ^^^^           ^^^^
            #           Must start with "data:"      Two newlines required!

            # Wait before next update
            await asyncio.sleep(0.5)

        # Send final event
        final_data = {"status": "completed", ...}
        yield f"data: {json.dumps(final_data)}\n\n"

    # Return streaming response
    return StreamingResponse(
        event_generator(),

```

```
        media_type="text/event-stream" # CRITICAL!
    )
```

## Frontend (JavaScript):

```
// Create EventSource connection
const eventSource = new EventSource('/api/upload/abc-123/stream');

// Handle incoming events
eventSource.onmessage = (event) => {
  const data = JSON.parse(event.data);
  console.log('Received:', data);

  // Update UI based on data
  updateUI(data);

  // Close when done
  if (data.status === 'completed') {
    eventSource.close();
  }
};

// Handle errors
eventSource.onerror = (error) => {
  console.error('SSE error:', error);
  eventSource.close();
};
```

## SSE Message Format

### HTTP Response:

```
HTTP/1.1 200 OK
Content-Type: text/event-stream
Cache-Control: no-cache
Connection: keep-alive

data: {"current_step":1,"step_status":"File uploaded"}\n\n
data: {"current_step":2,"step_status":"Parsing CSV..."}\n\n
data: {"current_step":3,"step_status":"Selected 10 patients"}\n\n
```

### Critical Requirements:

1. Content-Type MUST be **text/event-stream**
  2. Each event MUST start with **data:**
  3. Each event MUST end with **\n\n** (two newlines)
  4. Connection MUST stay open (keep-alive)
-

# Integration Guide

## Complete Frontend Integration Example

### Step 1: Upload CSV File

```
// HTML
<input type="file" id="csvFile" accept=".csv,.xlsx">
<select id="triggerEvent">
  <option value="ADT-A01">ADT-A01 - Register Patient</option>
  <option value="ADT-A04">ADT-A04 - Register Outpatient</option>
</select>
<label>
  <input type="checkbox" id="sendToMirth"> Send to Mirth
</label>
<button onclick="uploadFile()">Upload</button>

// JavaScript
async function uploadFile() {
  const fileInput = document.getElementById('csvFile');
  const triggerEvent = document.getElementById('triggerEvent').value;
  const sendToMirth = document.getElementById('sendToMirth').checked;

  if (!fileInput.files.length) {
    alert('Please select a file');
    return;
  }

  const formData = new FormData();
  formData.append('file', fileInput.files[0]);
  formData.append('trigger_event', triggerEvent);
  formData.append('send_to_mirth', sendToMirth);

  try {
    const response = await fetch('http://localhost:8000/api/upload', {
      method: 'POST',
      body: formData
    });

    if (!response.ok) {
      throw new Error(`HTTP error! status: ${response.status}`);
    }

    const data = await response.json();
    console.log('Upload started:', data.upload_id);

    // Start real-time progress tracking
    connectToProgressStream(data.upload_id);
  } catch (error) {
    console.error('Upload failed:', error);
    alert('Upload failed: ' + error.message);
  }
}
```

### Step 2: Connect to Progress Stream

```
// HTML
<div id="wizardSteps">
  <div class="step" data-step="1">1. Upload CSV</div>
  <div class="step" data-step="2">2. Parse Data</div>
```

```

<div class="step" data-step="3">3. Select Patients</div>
<div class="step" data-step="4">4. Create HL7</div>
<div class="step" data-step="5">5. Push to EMR</div>
<div class="step" data-step="6">6. Complete</div>
</div>
<div id="statusMessage"></div>
<div id="progressBar">
  <div id="progressFill" style="width: 0%></div>
  <span id="progressText">0/0</span>
</div>

// JavaScript
function connectToProgressStream(uploadId) {
  const eventSource = new EventSource(`http://localhost:8000/api/upload/${uploadId}/stream`);

  eventSource.onmessage = (event) => {
    const data = JSON.parse(event.data);
    console.log('Progress update:', data);

    // Update wizard steps
    highlightStep(data.current_step);

    // Update status message
    document.getElementById('statusMessage').textContent = data.step_status;

    // Update progress bar (if in step 4)
    if (data.generated_count && data.total_patients) {
      const percent = (data.generated_count / data.total_patients) * 100;
      document.getElementById('progressFill').style.width = percent + '%';
      document.getElementById('progressText').textContent =
        `${data.generated_count}/${data.total_patients}`;
    }

    // When complete, close connection and load results
    if (data.status === 'completed') {
      eventSource.close();
      console.log('Processing complete!');
      loadResults(uploadId);
    }
  }

  // Handle errors
  if (data.status === 'error') {
    eventSource.close();
    alert('Error: ' + data.error);
  }
};

eventSource.onerror = (error) => {
  console.error('SSE connection error:', error);
  eventSource.close();
  alert('Connection lost. Please refresh and try again.');
};

function highlightStep(stepNumber) {
  // Remove active class from all steps
  document.querySelectorAll('.step').forEach(el => {
    el.classList.remove('active');
  });

  // Add active class to current step
  const currentStep = document.querySelector(`[data-step="${stepNumber}"]`);
  if (currentStep) {
    currentStep.classList.add('active');
  }
}

```

### Step 3: Load Final Results

```

// HTML
<div id="resultsScreen" style="display:none">
  <h2>Processing Complete!</h2>
  <div class="summary">
    <div>Total Patients: <span id="totalPatients"></span></div>
    <div>Successful: <span id="successfulCount"></span></div>
    <div>Failed: <span id="failedCount"></span></div>
    <div>Sent to Mirth: <span id="mirthSentCount"></span></div>
  </div>
  <table id="resultsTable">
    <thead>
      <tr>
        <th>Patient Name</th>
        <th>Patient ID</th>
        <th>Status</th>
        <th>Sent to Mirth</th>
        <th>Actions</th>
      </tr>
    </thead>
    <tbody id="resultsTableBody"></tbody>
  </table>
</div>

// JavaScript
async function loadResults(uploadId) {
  try {
    const response = await fetch(`http://localhost:8000/api/upload/${uploadId}/results`);
    const results = await response.json();

    // Show results screen
    document.getElementById('resultsScreen').style.display = 'block';

    // Update summary
    document.getElementById('totalPatients').textContent = results.total_patients;
    document.getElementById('successfulCount').textContent = results.successful;
    document.getElementById('failedCount').textContent = results.failed;
    document.getElementById('mirthSentCount').textContent = results.mirth_successful;

    // Populate results table
    const tbody = document.getElementById('resultsTableBody');
    tbody.innerHTML = ''; // Clear existing rows

    results.messages.forEach(msg => {
      const row = tbody.insertRow();

      // Patient name
      const nameCell = row.insertCell();
      nameCell.textContent = msg.patient_name;

      // Patient ID
      const idCell = row.insertCell();
      idCell.textContent = msg.patient_id;

      // Status
      const statusCell = row.insertCell();
      statusCell.textContent = msg.status === 'success' ? '✓ Success' : '✗ Failed';
      statusCell.className = msg.status === 'success' ? 'success' : 'error';

      // Sent to Mirth
      const mirthCell = row.insertCell();
      mirthCell.textContent = msg.mirth_sent ? 'Yes' : 'No';

      // Actions
      const actionsCell = row.insertCell();
      const viewButton = document.createElement('button');
      viewButton.textContent = 'View HL7';
      viewButton.onclick = () => showHL7Message(msg.hl7_message);
      actionsCell.appendChild(viewButton);
    });
  } catch (error) {
    console.error('Failed to load results:', error);
    alert('Failed to load results: ' + error.message);
  }
}

```

```

        }

    }

function showHL7Message(hl7Message) {
    // Show HL7 message in modal or popup
    const modal = document.createElement('div');
    modal.className = 'modal';
    modal.innerHTML = `
        <div class="modal-content">
            <h3>HL7 Message</h3>
            <pre>${hl7Message}</pre>
            <button onclick="this.parentElement.parentElement.remove()">Close</button>
        </div>
    `;
    document.body.appendChild(modal);
}

```

## Step 4: Update Dashboard

```

// HTML
<div id="dashboard">
    <div class="tile">
        <h3>Total Processed</h3>
        <div id="totalProcessed" class="count">0</div>
    </div>
    <div class="tile">
        <h3>HL7 Messages</h3>
        <div id="hl7Messages" class="count">0</div>
    </div>
    <div class="tile">
        <h3>Success Rate</h3>
        <div id="successRate" class="count">0%</div>
    </div>
    <div class="system-status">
        <div>
            Mirth: <span id="mirthStatus" class="status-dot"></span>
        </div>
        <div>
            HL7 Parser: <span id="parserStatus" class="status-dot"></span>
        </div>
    </div>
</div>

// JavaScript
// Poll dashboard stats every 5 seconds
setInterval(async () => {
    try {
        const stats = await fetch('http://localhost:8000/api/dashboard/stats')
            .then(r => r.json());

        document.getElementById('totalProcessed').textContent = stats.total_processed;
        document.getElementById('hl7Messages').textContent = stats.hl7_messages;
        document.getElementById('successRate').textContent = stats.success_rate + '%';

    } catch (error) {
        console.error('Failed to fetch dashboard stats:', error);
    }
}, 5000);

// Poll system status every 10 seconds
setInterval(async () => {
    try {
        const status = await fetch('http://localhost:8000/api/dashboard/system-status')
            .then(r => r.json());

        // Update Mirth status indicator
        const mirthDot = document.getElementById('mirthStatus');
        if (status.openemr_connection.status === 'Active') {
            mirthDot.className = 'status-dot green';
        } else {
            mirthDot.className = 'status-dot red';
        }
    } catch (error) {
        console.error('Failed to fetch system status:', error);
    }
}, 10000);

```

```
        mirthDot.title = 'Mirth Connected';
    } else {
        mirthDot.className = 'status-dot red';
        mirthDot.title = 'Mirth Offline';
    }

    // Update parser status
    const parserDot = document.getElementById('parserStatus');
    if (status.hl7_parser.status === 'Running') {
        parserDot.className = 'status-dot green';
        parserDot.title = 'OpenAI Enabled';
    } else {
        parserDot.className = 'status-dot yellow';
        parserDot.title = 'Using Fallback Generator';
    }

} catch (error) {
    console.error('Failed to fetch system status:', error);
}
}, 10000);
```

---

## Testing Guide

### Testing in Swagger UI

#### Step 1: Start the Server

```
cd /Users/nagarajm/Work/SG/interface-wizard/actual-code
python main_with_fastapi.py --api
```

#### Step 2: Open Swagger UI

Navigate to: <http://localhost:8000/docs>

#### Step 3: Test Dashboard Stats

1. Find **Dashboard** section
2. Click on `GET /api/dashboard/stats`
3. Click "Try it out"
4. Click "Execute"
5. See response (should show all zeros initially)

## **Step 4: Test System Status**

1. Find `GET /api/dashboard/system-status`
2. Click "Try it out"
3. Click "Execute"
4. Check `openemr_connection.status` - should show "Active" if Mirth is running

## **Step 5: Test Real-Time Upload**

1. Find **Real-Time Upload** section
2. Click on `POST /api/upload`
3. Click "Try it out"
4. Click "Choose File" and select `test_patients.csv`
5. Set `trigger_event = "ADT-A01"`
6. Set `sendtomirth = true`
7. Click "Execute"
8. Copy the `upload_id` from response

## **Step 6: Test SSE Stream (Can't test in Swagger - use cURL)**

```
# Replace with your actual upload_id
curl -N "http://localhost:8000/api/upload/abc-123-def-456/stream"
```

You'll see real-time events stream in your terminal!

## **Step 7: Test Final Results**

1. Wait for processing to complete (~15 seconds for 5 patients)
2. Find `GET /api/upload/{upload_id}/results`
3. Click "Try it out"
4. Enter your `upload_id`
5. Click "Execute"
6. See complete results with all messages

---

## **Testing with cURL**

### **Test Dashboard Stats:**

```
curl http://localhost:8000/api/dashboard/stats
```

### **Test Upload:**

```
curl -X POST "http://localhost:8000/api/upload" \
-F "file=@test_patients.csv" \
-F "trigger_event=ADT-A01" \
-F "send_to_mirth=true"
```

### **Test SSE Stream:**

```
# Use -N flag for no buffering
curl -N "http://localhost:8000/api/upload/{upload_id}/stream"
```

### **Test Results:**

```
curl "http://localhost:8000/api/upload/{upload_id}/results"
```

---

## Troubleshooting

### Common Issues

#### 1. "Connection refused" on Mirth

##### **Symptom:**

Connection refused! Check Mirth is running and channel is listening.

##### **Causes:**

- Mirth Connect not running
- Channel not deployed/started
- Wrong port (should be 6661)

**Solution:**

1. Open Mirth Connect Administrator (<https://localhost:8443>)
2. Go to Channels tab
3. Verify "IW Listener HL7" channel exists
4. Check status is "Started" (green)
5. Verify listener port is 6661

**Verify Connection:**

```
# Test if port 6661 is listening  
netstat -an | grep 6661  
  
# Or telnet test  
telnet localhost 6661
```

---

**2. "Upload session not found"****Symptom:**

```
{  
    "detail": "Upload session not found"  
}
```

**Causes:**

- Invalid upload\_id
- Server restarted (sessions lost)
- Session expired (unlikely with in-memory)

**Solution:**

1. Verify you copied the correct upload\_id
  2. Check if server restarted (sessions are in-memory)
  3. Upload file again to get new upload\_id
- 

**3. Dashboard stats show 0**

**Symptom:**

All dashboard stats are 0

**Causes:**

- No uploads processed yet
- Server restarted (stats reset)

**Solution:**

This is normal! Stats will update after processing uploads.

---

## 4. SSE stream not connecting

**Symptom:**

EventSource connection fails or times out

**Causes:**

- CORS issues
- Invalid upload\_id
- Browser doesn't support SSE (rare)

**Solution:****Check CORS:**

```
# In main_with_fastapi.py, add CORS middleware if needed
from fastapi.middleware.cors import CORSMiddleware

app.add_middleware(
    CORSMiddleware,
    allow_origins=["*"], # Or specific origins
    allow_credentials=True,
    allow_methods=["*"],
    allow_headers=["*"],
)
```

**Use polling instead:**

```
// If SSE doesn't work, use polling
```

```
async function pollStatus(uploadId) {
  const interval = setInterval(async () => {
    const status = await fetch(`/api/upload/${uploadId}/status`)
      .then(r => r.json());

    if (status.status === 'completed') {
      clearInterval(interval);
    }
  }, 1000);
}
```

---

## 5. OpenAI API errors

### Symptom:

■ Remote generation failed: ... Falling back to local generator.

### Causes:

- No API key configured
- API key invalid
- Rate limit exceeded
- Network error

### Solution:

#### Check API key:

```
# In main_with_fastapi.py, line 46
OPENAI_API_KEY = "sk-..." # Replace with valid key
```

#### Fallback works fine:

The system automatically uses the deterministic fallback generator if OpenAI fails. Messages will still be valid HL7 v2.5 format!

---

## Glossary

## **Healthcare/HL7 Terms**

### **HL7 (Health Level 7)**

- International standard for healthcare data exchange
- Version 2.x is pipe-delimited text format
- Version 3.x and FHIR are XML/JSON based

### **ADT (Admission/Discharge/Transfer)**

- Message type for patient movement
- ADT-A01 = Admit/Register Patient
- ADT-A04 = Register Outpatient
- ADT-A08 = Update Patient Information

### **MLLP (Minimal Lower Layer Protocol)**

- Protocol for sending HL7 messages over TCP/IP
- Wraps message with special bytes:
- Start: (0x0B)
- End: (0x1C 0x0D)

### **Mirth Connect**

- Open-source HL7 message router
- Receives, transforms, and routes HL7 messages
- Has channels with source/destination connectors

### **OpenEMR**

- Open-source Electronic Medical Record system
- Uses MySQL database
- Stores patient data in `patient_data` table

### **ACK (Acknowledgment)**

- Response message confirming receipt
- AA = Application Accept (success)
- AE = Application Error (rejected)

## Technical Terms

### FastAPI

- Modern Python web framework
- Built on Starlette and Pydantic
- Automatic API documentation (Swagger/OpenAPI)
- Async/await support

### Pydantic

- Data validation using Python type hints
- Automatic JSON serialization
- Runtime type checking

### Server-Sent Events (SSE)

- Standard for server-to-client real-time updates
- Uses HTTP (not WebSocket)
- One-way communication
- Auto-reconnect

### Async/Await

- Python concurrency feature
- Non-blocking I/O operations
- Allows multiple tasks to run concurrently

### UUID (Universally Unique Identifier)

- 128-bit unique identifier
- Example: `abc-123-def-456`
- Used for session IDs

### In-Memory Storage

- Data stored in RAM (Python dictionaries)
- Fast access
- Lost on restart
- Alternative: Database persistence

## Rate Limiting

- Intentional delay between operations
- Prevents overwhelming services
- Example: 1 second between HL7 generations

## CORS (Cross-Origin Resource Sharing)

- Security feature for web browsers
- Controls which domains can access API
- Needed for frontend-backend communication

---

## Appendix: Example HL7 Message

```
MSH|^~\&|SMART_APP|SMART_FAC|REC_APP|REC_FAC|20251220140030||ADT^A01|abc-123-def-456|P|2.5
EVN|A01|20251220140030
PID|1||PAT001||Doe^John||19900115|M||123 Main St, Boston, MA 02101
PV1|1|I
```

### Segment Breakdown:

#### MSH (Message Header):

- Field 1: **MSH** - Segment type
- Field 2: **^~\&** - Encoding characters
- Field 3: **SMART\_APP** - Sending application
- Field 4: **SMART\_FAC** - Sending facility
- Field 5: **REC\_APP** - Receiving application
- Field 6: **REC\_FAC** - Receiving facility
- Field 7: **20251220140030** - Timestamp (YYYYMMDDHHmmss)
- Field 9: **ADT^A01** - Message type (ADT) and trigger event (A01)
- Field 10: **abc-123-def-456** - Message control ID (unique)
- Field 11: **P** - Processing ID (P = Production)
- Field 12: **2.5** - HL7 version

#### EVN (Event):

- Field 1: `A01` - Event type code
- Field 2: `20251220140030` - Event timestamp

#### PID (Patient Identification):

- Field 1: `1` - Set ID
- Field 3: `PAT001` - Patient identifier (MRN)
- Field 5: `Doe^John` - Patient name (Last^First)
- Field 7: `19900115` - Date of birth (YYYYMMDD)
- Field 8: `M` - Gender (M/F/U)
- Field 11: `123 Main St, Boston, MA 02101` - Address

#### PV1 (Patient Visit):

- Field 1: `1` - Set ID
  - Field 2: `I` - Patient class (I = Inpatient, O = Outpatient)
- 

## Conclusion

This documentation provides a complete guide to the Interface Wizard Backend API. Key takeaways:

- **6 new endpoints** for real-time CSV processing and dashboard statistics
- **Server-Sent Events (SSE)** for live progress updates
- **In-memory storage** for fast, simple session management
- **Mirth Connect integration** via MLLP protocol
- **OpenAI-powered HL7 generation** with fallback
- **Complete Swagger documentation** at `/docs`
- **Production-ready** with proper error handling and validation

## Next Steps

**For Backend Developers:**

1. Review code structure in `mainwithfastapi.py`
2. Test endpoints in Swagger UI
3. Understand SSE implementation
4. Customize as needed for your environment

**For Frontend Developers:**

1. Review Integration Guide section
2. Test API endpoints from your framework
3. Implement SSE connection
4. Build 6-step wizard UI

**For System Administrators:**

1. Deploy FastAPI application
2. Configure Mirth Connect
3. Set up OpenAI API key
4. Monitor logs and performance

**For QA Engineers:**

1. Follow Testing Guide
2. Test with sample CSV files
3. Verify Mirth integration
4. Test error scenarios

---

**Document Version:** 1.0

**Last Updated:** December 20, 2025

**For Questions:** Contact Interface Wizard Development Team